

R-60S Analyzer



R-60S is an advanced semi-automatic chemistry analyzer for small and medium clinics and laboratories.

This compact urine analyzer is well designed for easy use and a user-friendly system with a built in printer.

USB barcode scanner is optional

BioSystems Diagnostics Private Limited

A4, Sipcot Industrial Park | Irungattukottai | TN-602105

+91-44-67121911 | support@biosystems.in | www.biosystems.in

R-60S Analyzer



Technical features

- Reflectance photometer for urinalysis
- 50 samples/hour
- Memory : 2,000 pcs
- 4.3 inch touch LCD
- Built-in thermal printer(57 x 50mm paper)
- Compact size (235 x 2232 x 85 mm)
- Lightweight, only 840g
- Input : 100-240V/Output : DC12V 3.33A
- Interface : USB type B
- Optional Accessory : USB Barcode Scanner
- Auto- Calibration and System Check control
- Read up to 11 parameters



Ordering Information

Brand	Parameter	Cat.No
CYBOW R-60S	CYBOW 2AC, 10, 11	5A00
Combostik R-60S	Combostik 2AC, 10, 11	5A40
DUS R-60S	DUS 2AC, 10, 11	5A60
UriDoctor Vet Plus	Vet 10, 11C, 11AC, BHB, LDH, MUN	5B90
UriDoctor Plus	UriDoctor 2AC, 11	5A80

BioSystems Diagnostics Private Limited

A4, Sipcot Industrial Park | Irungattukottai | TN-602105

+91-44-67121911 | support@biosystems.in | www.biosystems.in