



Product
List 2026

Clinical analysis

BioSystems Sujetsterms

Tests sorted alphabetically

Test	Code	RA	RB	mL/Kit	Comment
α1-ACID GLYCOPROTEIN					
<i>Immunoturbidimetric assay</i>					
Manual	31928	1 x 50 mL	-	50	
BA Automated System	23107	1 x 60 mL	1 x 15 mL	75	
α-AMYLASE DIRECT					
<i>Direct Substrate</i>					
Manual	11583	5 x 5 mL	-	25	
	11550	6 x 25 mL	-	150	
A15/A25 Automated System	12550	5 x 20 mL	-	100	
BA Automated System	21550	8 x 20 mL	-	160	
	23550	4 x 20 mL	-	80	
α-AMYLASE-EPS					
<i>IFCC</i>					
Manual	11534	1 x 32 mL	1 x 8 mL	40	WR
A15/A25 Automated System	12535	1 x 32 mL	1 x 8 mL	40	WR
BA Automated System	21534	2 x 60 mL	2 x 15 mL	150	
α-AMYLASE PANCREATIC					
<i>Immuno-inhibition</i>					
Manual	11799	1 x 20 mL	1 x 5 mL	25	
A15/A25 Automated System	12799	1 x 40 mL	1 x 10 mL	50	
BA Automated System	21799	2 x 60 mL	2 x 15 mL	150	
α-GLUCOSIDASE					
<i>Glucose oxidase/Peroxidase</i>					
A15/A25 Automated System	12522	2 x 20 mL	2 x 5 mL	50	Std Inc
BA Automated System	21522	2 x 20 mL	2 x 5 mL	50	Std Inc
ACID PHOSPHATASE (ACP)					
<i>Naphtyl phosphate/Pentanediol</i>					
Manual (Total + Non prostatic)	11548	4x10 mL	lyophilised	40	WR
ADENOSINE DEAMINASE (ADA)					
<i>Adenosine-Glutamate Dehydrogenase</i>					
Manual / A15/A25 Automated System	12754	4 x 8 mL	1 x 10 mL	40	WR
BA Automated System	23754	2 x 16 mL	1 x 10 mL	40	WR
ALANINE AMINOTRANSFERASE (ALT/GPT)					
<i>IFCC</i>					
Manual	11832	1 x 40 mL	1 x 10 mL	50	WR
	11533	1 x 160 mL	1 x 40 mL	200	WR
	11568	1 x 400 mL	1 x 100 mL	500	WR
	11562	1 x 800 mL	1 x 200 mL	1000	WR
A15/A25 Automated System	12533	5 x 40 mL	5 x 10 mL	250	WR
BA Automated System	21533	8 x 60 mL	8 x 15 mL	600	
	23533	4 x 60 mL	4 x 15 mL	300	
Pyridoxal Phosphate (ALT)	11667	1 x 5 mL			
ALBUMIN					
<i>Bromocresol Green</i>					
Manual	11573	1 x 250 mL	-	250	Std Inc
	11547	2 x 250 mL	-	500	Std Inc
A15/A25 Automated System	12547	5 x 50 mL	-	250	
BA Automated System	21547	10 x 60 mL	-	600	
	23547	4 x 60 mL	-	240	

In bold italics: methods | **WR**: Working Reagent | **RA**: Reagent A Preparation | **RB**: Reagent B Preparation | **Std Inc**: Standard Included

Test	Code	RA	RB	mL/Kit	Comment
ALBUMIN (MICROALBUMINURIA)					
<i>Immunoturbidimetric latex assay</i>					
Manual	31324	1 x 16 mL	1 x 4 mL	20	WR
	31924	1 x 40 mL	1 x 10 mL	50	WR
A15/A25 Automated System	13324	1 x 40 mL	1 x 10 mL	50	WR
BA Automated System	22324	4 x 60 mL	4 x 15 mL	300	
	23324	1 x 60 mL	1 x 15 mL	75	
ALKALINE PHOSPHATASE (ALP)-AMP					
<i>AMP Buffer (IFCC)</i>					
Manual	11592	1 x 40 mL	1 x 10 mL	50	WR
	11593	1 x 160 mL	1 x 40 mL	200	WR
	11598	4 x 100 mL	2 x 50 mL	500	WR
A15/A25 Automated System	12518	5 x 16 mL	2 x 10 mL	100	WR
BA Automated System	21592	4 x 60 mL	4 x 15 mL	300	
	23592	1 x 60 mL	1 x 15 mL	75	
ALKALINE PHOSPHATASE (ALP)-DEA					
<i>Diethanolamine Buffer</i>					
Manual	11590	1 x 40 mL	1 x 10 mL	50	WR
	11591	4 x 40 mL	4 x 10 mL	200	WR
	11597	4 x 100 mL	2 x 50 mL	500	WR
A15/A25 Automated System	12514	5 x 16 mL	2 x 10 mL	100	WR
BA Automated System	21590	4 x 60 mL	4 x 15 mL	300	
	23590	1 x 60 mL	1 x 15 mL	75	
AMMONIA					
<i>Glutamate Dehydrogenase</i>					
Manual / A15/A25 Automated System	12532	1 x 20 mL	1 x 7 mL	27	RA
BA Automated System	23532	1 x 20 mL	1 x 7 mL	27	RA
ANGIOTENSIN CONVERTING ENZYME (ACE)					
<i>FAPGG</i>					
Manual / A15/A25 Automated System	12796	1 x 50 mL	-	50	
BA Automated System	21796	2 x 60 mL	-	120	
ANTI-STREPTOLYSIN O (ASO)					
<i>Immunoturbidimetric latex assay</i>					
Manual	31323	1 x 16 mL	1 x 4 mL	20	WR
	31923	1 x 40 mL	1 x 10 mL	50	WR
	31031	2 x 160 mL	2 x 40 mL	400	WR
A15/A25 Automated System	13923	1 x 40 mL	1 x 10 mL	50	WR
BA Automated System	22923	2 x 60 mL	2 x 15 mL	150	
	23923	1 x 60 mL	1 x 15 mL	75	
APOLIPOPROTEIN A-1 (APO A-1)					
<i>Immunoturbidimetric assay</i>					
Manual	31095	1 x 40 mL	1 x 10 mL	50	
BA Automated System	23095	1 x 60 mL	1 x 15 mL	75	
APOLIPOPROTEIN B (APO B)					
<i>Immunoturbidimetric assay</i>					
Manual	31098	1 x 40 mL	1 x 10 mL	50	
BA Automated System	23098	1 x 60 mL	1 x 15 mL	75	

Test	Code	RA	RB	mL/Kit	Comment
ASPARTATE AMINOTRANSFERASE (AST/GOT)					
<i>IFCC</i>					
Manual	11830	1 x 40 mL	1 x 10 mL	50	WR
	11531	1 x 160 mL	1 x 40 mL	200	WR
	11567	1 x 400 mL	1 x 100 mL	500	WR
	11561	1 x 800 mL	1 x 200 mL	1000	WR
A15/A25 Automated System	12531	5 x 40 mL	5 x 10 mL	250	WR
BA Automated System	21531	8 x 60 mL	8 x 15 mL	600	
	23531	4 x 60 mL	4 x 15 mL	300	
Pyridoxal Phosphate (AST)	11666	1 x 5 mL			
β-HYDROXYBUTYRATE					
<i>Hydroxybutyrate Dehydrogenase/Diaphorase</i>					
Manual/A15/A25 Automated System	12525	1 x 40 mL	1 x 10 mL	50	
BA Automated System	21525	1 x 60 mL	1 x 15 mL	75	
β2-MICROGLOBULIN					
<i>Immunoturbidimetric latex assay</i>					
Manual	31920	1 x 800 mL	1 x 200 mL	1000	WR
	31925	1 x 40 mL	1 x 10 mL	50	WR
BA Automated System	23925	1 x 60 mL	1 x 15 mL	75	
BILIRUBIN (DIRECT)					
<i>Diazotized Sulfanilic</i>					
Manual	11511	4 x 40 mL	4 x 10 mL	200	WR
	11545	2 x 400 mL	2 x 100 mL	1000	WR
BILIRUBIN (DIRECT) DPD					
<i>Dichlorophenyl Diazonium</i>					
A15/A25 Automated System	12504	5 x 40 mL	5 x 10 mL	250	
BA Automated System	21504	4 x 60 mL	4 x 15 mL	300	
	23504	1 x 60 mL	1 x 15 mL	75	
BILIRUBIN (TOTAL)					
<i>Diazotized Sulfanilic</i>					
Manual	11510	4 x 40 mL	4 x 10 mL	200	WR
	11544	2 x 400 mL	2 x 100 mL	1000	WR
BILIRUBIN (TOTAL) DPD					
<i>Dichlorophenyl Diazonium</i>					
A15/A25 Automated System	12506	5 x 40 mL	5 x 10 mL	250	
BA Automated System	21506	8 x 60 mL	8 x 15 mL	600	
	23506	4 x 60 mL	4 x 15 mL	300	
BILIRUBIN (TOTAL AND DIRECT)					
<i>Diazotized Sulfanilic</i>					
Manual	11515	2 x 50 mL	2 x 50 mL	200	WR
	11555	1 x 500 mL	1 x 500 mL	1000	WR
CALCIUM-ARSENZO					
<i>Arsenazo III</i>					
Manual	11570	1 x 200 mL	-	200	Std Inc
	11571	1 x 500 mL	-	500	Std Inc
A15/A25 Automated System	12570	10 x 50 mL	-	500	
BA Automated System	21570	10 x 60 mL	-	600	
	23570	4 x 60 mL	-	240	

Test	Code	RA	RB	mL/Kit	Comment
CALCIUM-CRESOLPHTHALEIN					
<i>O-Cresolphthalein Complexone</i>					
Manual	11811	1 x 160 mL	1 x 40 mL	200	Std Inc, WR
	11812	1 x 400 mL	2 x 50 mL	500	Std Inc, WR
A15/A25 Automated System	12513	5 x 40 mL	5 x 10 mL	250	
BA Automated System	21511	8 x 60 mL	8 x 15 mL	600	
	23511	2 x 60 mL	2 x 15 mL	150	
CALCIUM-MTB					
<i>Metihylthymol Blue</i>					
Manual	11527	2 x 50 mL	2 x 50 mL	200	Std Inc, WR
	11507	1 x 250 mL	1 x 250 mL	500	Std Inc, WR
CARBON DIOXIDE (CO₂)					
<i>Phosphoenolpyruvate Carboxylase/Malate Dehydrogenase</i>					
Manual/A15/A25 Automated System	11558	1 x 50 mL	-	50	
BA Automated System	21558	2 x 60 mL	-	120	
CERULOPLASMIN					
<i>Turbidimetry</i>					
A15/A25 Automated System	13340	1 x 40mL	1 x 10mL	50	
BA Automated System	22340	1 x 40mL	1 x 10mL	50	
CHOLESTEROL					
<i>Cholesterol Oxidase/Peroxidase</i>					
Manual	11805	1 x 50 mL	-	50	Std Inc
	11505	1 x 200 mL	-	200	Std Inc
	11506	1 x 500 mL	-	500	Std Inc
	11539	1 x 1000 mL	-	1000	Std Inc
A15/A25 Automated System	12505	10 x 50 mL	-	500	
BA Automated System	21505	10 x 60 mL	-	600	
	23505	4 x 60 mL	-	240	
CHOLESTEROL HDL					
<i>Phosphotungstate/ Mg²⁺-Cholesterol Oxidase/Peroxidase</i>					
Manual	11523	2 x 50 mL	2 x 50 mL	200	Std Inc
CHOLESTEROL HDL DIRECT					
<i>TOOS</i>					
Manual	11757	1x60mL	1x20mL	80	
A15/A25 Automated System	12757	3x20mL	1x20mL	80	
BA Automated System	21594	4x60mL	4x20mL	320	
	23594	1x60mL	1x20mL	80	
CHOLESTEROL HDL PRECIPITATING REAGENT					
<i>Phosphotungstate/ Mg²⁺</i>					
Manual	11648	1 x 50 mL	-	50	Std Inc
CHOLESTEROL LDL DIRECT					
<i>TOOS</i>					
Manual	11785	1x60mL	1x20mL	80	
A15/A25 Automated System	12785	3x20mL	1x20mL	80	
BA Automated System	21785	2x60mL	2x20mL	160	
	23785	1x60mL	1x20mL	80	
CHOLESTEROL LDL PRECIPITATING REAGENT					
<i>Polyvinyl Sulphate/ Polyethyleneglycol</i>					
Manual	11579	1 x 20 mL	-	20	

Test	Code	RA	RB	mL/Kit	Comment
CHOLINESTERASE (CHE)					
<i>Butyrylthiocholine</i>					
Manual	11588	1 x 40 mL	1 x 10 mL	50	WR
	11589	4 x 40 mL	4 x 10 mL	200	WR
A15/A25 Automated System	12588	1 x 40 mL	1 x 10 mL	50	WR
BA Automated System	21588	2 x 60 mL	2 x 15 mL	150	
CITRATE					
<i>Citrate Lyase/ Malate Dehydrogenase</i>					
Manual / A15/A25 Automated System	11895	1 x 20 mL	1 x 5 mL	25	RB + Std Inc
	11795	2 x 20 mL	2 x 5 mL	50	RB + Std Inc
BA Automated System	23795	2 x 20 mL	2 x 5 mL	50	RB + Std Inc
COMPLEMENT COMPONENT C3					
<i>Immunoturbidimetric assay</i>					
Manual	31084	1 x 20 mL	-	20	
	31073	1 x 50 mL	-	50	
	31079	1 x 250 mL	-	250	
A15/A25 Automated System	13084	1 x 50 mL	-	50	
BA Automated System	23103	1 x 60 mL	1 x 15 mL	75	
COMPLEMENT COMPONENT C4					
<i>Immunoturbidimetric assay</i>					
Manual	31085	1 x 20 mL	-	20	
	31074	1 x 50 mL	-	50	
	31080	1 x 250 mL	-	250	
A15/A25 Automated System	13085	1 x 50 mL	-	50	
BA Automated System	23104	1 x 60 mL	1 x 15 mL	75	
COPPER - PAESA					
<i>3,5-DiBr-PAESA</i>					
Manual	11837	1 x 40 mL	1 x 10 mL	50	RA + Std Inc
C-REACTIVE PROTEIN (CRP)					
<i>Immunoturbidimetric latex assay</i>					
Manual	31321	1 x 16 mL	1 x 4 mL	20	WR
	31921	1 x 40 mL	1 x 10 mL	50	WR
	31029	2 x 160 mL	2 x 40 mL	400	WR
A15/A25 Automated System	13921	2 x 40 mL	2 x 10 mL	100	WR
BA Automated System	22921	4 x 60 mL	4 x 15 mL	300	
	23921	1 x 60 mL	1 x 15 mL	75	
C-REACTIVE PROTEIN hs (CRP-hs)					
<i>Immunoturbidimetric latex assay, High Sensitivity</i>					
Manual	31927	1 x 40 mL	1 x 10 mL	50	WR
A15/A25 Automated System	13927	1 x 40 mL	1 x 10 mL	50	WR
BA Automated System	22927	2 x 60 mL	2 x 15 mL	150	
	23927	1 x 60 mL	1 x 15 mL	75	
C-REACTIVE PROTEIN-er (CRP-er)					
<i>Immunoturbidimetric latex assay, Extended Range</i>					
BA Automated System	22942	2 x 60 mL	2 x 15 mL	150	
	23942	1 x 60 mL	1 x 15 mL	75	

Test	Code	RA	RB	mL/Kit	Comment
CREATINE KINASE (CK)					
<i>IFCC</i>					
Manual	11790	1 x 40 mL	1 x 10 mL	50	WR
	11791	4 x 40 mL	4 x 10 mL	200	WR
A15/A25 Automated System	12524	3 x 12 mL	1 x 10 mL	45	WR
BA Automated System	21790	2 x 60 mL	2 x 15 mL	150	
	23790	1 x 60 mL	1 x 15 mL	75	
CREATINE KINASE-MB (CK-MB)					
<i>Immunoinhibition</i>					
Manual	11792	1 x 40 mL	1 x 10 mL	50	WR
A15/A25 Automated System	12566	3 x 12 mL	1 x 10 mL	45	WR
BA Automated System	21792	2 x 60 mL	2 x 15 mL	150	
	23792	1 x 60 mL	1 x 15 mL	75	
CREATININE					
<i>Jaffé-compensated</i>					
Manual	11802	1 x 50 mL	1 x 50 mL	100	Std Inc, WR
	11502	2 x 50 mL	2 x 50 mL	200	Std Inc, WR
	11542	1 x 500 mL	1 x 500 mL	1000	Std Inc, WR
A15/A25 Automated System	12502	5 x 50 mL	5 x 50 mL	500	WR
BA Automated System	21502	5 x 60 mL	5 x 60 mL	600	
	23502	5 x 20 mL	5 x 20 mL	200	
CREATININE ENZYMATIC					
<i>Enzymatic</i>					
Manual	11734	1 x 60 mL	1 x 20 mL	80	Std Inc
A15/A25 Automated System	12734	1 x 45 mL	1 x 15 mL	60	
BA Automated System	21734	2 x 60 mL	2 x 20 mL	160	
D-DIMER (DD)					
<i>Immunoturbidimetric latex assay</i>					
Manual / A15/A25 Automated System	13519	1 x 20 mL	1 x 10 mL	30	
BA Automated System	21519	1 x 20 mL	1 x 10 mL	30	
ETHANOL					
<i>Alcohol Dehydrogenase</i>					
Manual / A15/A25 Automated System	12789	1 x 20 mL	1 x 7 mL	27	
BA Automated System	21789	2 x 20 mL	2 x 7 mL	54	
FECAL CALPROTECTIN (FCAL)					
<i>Turbidimetry</i>					
A15S Automated System	12330	1 x 40 mL	1 x 10 mL	50	Reagent C included (1 x 10 mL)
STOOL Collection Tube	31329	-	-	50 tubes	
Open STOOL Collection Tube (oSCT)	31334	-	-	50 tubes	
FECAL OCCULT BLOOD (FOB)					
<i>Turbidimetry</i>					
A15S Automated System	12325	1 x 40 mL	1 x 10 mL	50	Reagent C included (1 x 10 mL)
STOOL Collection Tube	31329	-	-	50 tubes	
Open STOOL Collection Tube (oSCT)	31334	-	-	50 tubes	

Test	Code	RA	RB	mL/Kit	Comment
FERRITIN					
<i>Immunoturbidimetric latex assay</i>					
Manual	31934	1 x 10 mL	1 x 5 mL	15	WR
	31935	1 x 30 mL	1 x 15 mL	45	WR
A15/A25 Automated System	13934	1 x 30 mL	1 x 15 mL	45	WR
BA Automated System	22934	2 x 40 mL	2 x 20 mL	120	
	23934	1 x 40 mL	1 x 20 mL	60	
FIBRINOGEN					
<i>Immunoturbidimetric assay</i>					
Manual	31600	1 x 40 mL	1 x 10 mL	50	
A15/A25 Automated System	13600	1 x 40 mL	1 x 10 mL	50	
BA Automated System	23804	1 x 60 mL	1 x 15 mL	75	
FRUCTOSAMINE					
<i>NBT</i>					
Manual	11046	2 x 50 mL	-	100	Std Inc
FRUCTOSE					
<i>Hexokinase/ Phosphoglucose Isomerase</i>					
Manual / A15/A25 Automated System	11794	1 x 40 mL	1 x 10 mL	50	Std Inc, WR
BA Automated System	23794	1 x 60 mL	1 x 15 mL	75	Std Inc
γ-GLUTAMYL TRANSFERASE (γ-GT)					
<i>IFCC</i>					
Manual	11584	1 x 40 mL	1 x 10 mL	50	WR
	11520	1 x 160 mL	1 x 40 mL	200	WR
A15/A25 Automated System	12520	5 x 40 mL	5 x 10 mL	250	WR
BA Automated System	21520	4 x 60 mL	4 x 15 mL	300	
	23520	2 x 60 mL	2 x 15 mL	150	
GLUCOSE					
<i>Glucose Oxidase/ Peroxidase</i>					
Manual	11803	1 x 50 mL	-	50	Std Inc
	11503	1 x 200 mL	-	200	Std Inc
	11504	1 x 500 mL	-	500	Std Inc
	11538	1 x 1000 mL	-	1000	Std Inc
A15/A25 Automated System	12503	10 x 50 mL	-	500	
BA Automated System	21503	10 x 60 mL	-	600	
	23503	4 x 60 mL	-	240	
GLUCOSE-6-PHOSPHATE DEHYDROGENASE					
<i>Glucose-6-Phosphate</i>					
Manual / A15/A25 Automated System	12603	1 x 40 mL	1 x 10 mL	50	Std Inc
BA Automated System	21603	1 x 60 mL	1 x 15 mL	75	Std Inc
GLUCOSE-HEXOKINASE					
<i>Hexokinase</i>					
Manual	11656	1 x 160 mL	1 x 40 mL	200	Std Inc, WR
A15/A25 Automated System	12756	2 x 40 mL	2 x 10 mL	100	
BA Automated System	21656	4 x 60 mL	4 x 15 mL	300	
	23656	1 x 60 mL	1 x 15 mL	75	
HAPTOGLOBIN					
<i>Immunoturbidimetric assay</i>					
Manual / A15/A25 Automated System	13218	1 x 40 mL	1 x 10 mL	50	
BA Automated System	23218	1 x 40 mL	1 x 10 mL	50	

Test	Code	RA	RB	mL/Kit	Comment
HEMOGLOBIN A1C-DIRECT (HbA1C-DIR)					
<i>Immunoturbidimetric latex assay</i>					
Manual	31047	1 x 50 mL	1 x 10 mL	60	
A15/A25 Automated System	13047	1 x 50 mL	1 x 10 mL	60	
BA Automated System	22047	2 x 60 mL	2 x 12 mL	144	
	22147	1 x 60 mL	1 x 12 mL	72	
HOMOCYSTEINE					
<i>Enzymatic Cycling</i>					
A15/A25 Automated System	12737	1 x 40 mL	1 x 10.8 mL	50.8	
BA Automated System	23737	1 x 20 mL	1 x 5.4 mL	25.4	
IMMUNOGLOBULIN A (IgA)					
<i>Immunoturbidimetric assay</i>					
Manual	31082	1 x 20 mL	-	20	
	31071	1 x 50 mL	-	50	
	31077	1 x 250 mL	-	250	
A15/A25 Automated System	13082	1 x 50 mL	-	50	
BA Automated System	23101	1 x 60 mL	1 x 15 mL	75	
IMMUNOGLOBULIN G (IgG)					
<i>Immunoturbidimetric assay</i>					
Manual	31081	1 x 20 mL	-	20	
	31070	1 x 50 mL	-	50	
	31076	1 x 250 mL	-	250	
A15/A25 Automated System	13081	1 x 50 mL	-	50	
BA Automated System	23100	1 x 60 mL	1 x 15 mL	75	
IMMUNOGLOBULIN M (IgM)					
<i>Immunoturbidimetric assay</i>					
Manual	31083	1 x 20 mL	-	20	
	31072	1 x 50 mL	-	50	
	31078	1 x 250 mL	-	250	
A15/A25 Automated System	13083	1 x 50 mL	-	50	
BA Automated System	23102	1 x 60 mL	1 x 15 mL	75	
IRON-CHROMAZUROL					
<i>Chromazurol B</i>					
Manual	11546	4 x 50 mL	-	200	Std Inc
IRON-FERROZINE					
<i>Ferrozine</i>					
Manual	11509	4 x 40 mL	4 x 10 mL	200	Std Inc, WR
A15/A25 Automated System	12509	5 x 40 mL	5 x 10 mL	250	
BA Automated System	21509	4 x 60 mL	4 x 15 mL	300	
	23509	2 x 60 mL	2 x 15 mL	150	
LACTATE					
<i>Lactate Oxidase/ Peroxidase</i>					
Manual	11736	1 x 40 mL	1 x 10 mL	50	WR
A15/A25 Automated System	12736	2 x 40 mL	2 x 10 mL	100	
BA Automated System	21736	4 x 60 mL	4 x 15 mL	300	
	23736	1 x 60 mL	1 x 15 mL	75	

Test	Code	RA	RB	mL/Kit	Comment
LACTATE DEHYDROGENASE (LDH)					
<i>Pyruvate</i>					
Manual	11580	1 x 40 mL	1 x 10 mL	50	WR
	11581	1 x 160 mL	1 x 40 mL	200	WR
A15/A25 Automated System	12580	5 x 40 mL	5 x 10 mL	250	WR
BA Automated System	21580	8 x 60 mL	8 x 15 mL	600	
	23580	4 x 60 mL	4 x 15 mL	300	
LACTATE DEHYDROGENASE (LDH)-IFCC					
<i>IFCC</i>					
Manual/ A15/A25 Automated System	11587	1 x 160 mL	1 x 40 mL	200	WR (Manual)
	11586	1 x 40 mL	1 x 10 mL	50	WR
BA Automated System	21586	8 x 60 mL	8 x 15 mL	600	
	23586	2 x 60 mL	2 x 15 mL	150	
LIPASE					
<i>DGGR</i>					
Manual	11760	1 x 20 mL	1 x 10 mL	30	Std Inc
A15/A25 Automated System	12760	1 x 20 mL	1 x 10 mL	30	
BA Automated System	21760	1 x 20 mL	1 x 10 mL	30	
MAGNESIUM					
<i>Xylydyl Blue</i>					
Manual	11797	4 x 40 mL	4 x 10 mL	200	Std Inc, WR
A15/A25 Automated System	12797	5 x 16 mL	2 x 10 mL	100	WR
BA Automated System	21797	2 x 60 mL	2 x 15 mL	150	
	23797	1 x 60 mL	1 x 15 mL	75	
NON-ESTERIFIED FATTY ACIDS (NEFA)					
<i>Acyl-CoA Oxidase/Peroxidase</i>					
A15/A25 Automated System	12540	2 x 20 mL	2 x 7 mL	54	RA, RB
BA Automated System	21840	2 x 20 mL	2 x 7 mL	54	RA, RB
OXALATE					
<i>Oxalate Oxidase/Peroxidase</i>					
Manual / A15/A25 Automated System	12539	1 x 20 mL	1 x 5 mL	1 x 25 mL R1 + 20 tubes (pretreatment reagent)	Std Inc, RB
Pretreatment Reagents	11839	-	-	5 x 25 mL R1 + 100 tubes	
BA Automated System (w/o tubes and pretreatment reagent)	23539	1 x 20 mL	1 x 5 mL	25	Std Inc, RB
PHOSPHORUS					
<i>Phosphomolybdate/UV</i>					
Manual	11508	3 x 40 mL	1 x 50 mL	170	Std Inc, WR
A15/A25 Automated System	12508	3 x 24 mL	2 x 15 mL	100	
BA Automated System	21518	4 x 50 mL	4 x 20 mL	280	
	23518	1 x 50 mL	1 x 20 mL	70	
PREALBUMIN					
<i>Immunoturbidimetric assay</i>					
Manual	31929	1 x 50 mL	-	50	
BA Automated System	23106	1 x 60 mL	1 x 15 mL	75	

Test	Code	RA	RB	mL/Kit	Comment
PROTEIN (TOTAL)					
<i>Biuret</i>					
Manual	11800	1 x 50 mL	-	50	Std Inc
	11572	1 x 250 mL	-	250	Std Inc
	11500	2 x 250 mL	-	500	Std Inc
	11553	1 x 1000 mL	-	1000	Std Inc
A15/A25 Automated System	12500	10 x 50 mL	-	500	
BA Automated System	21513	2 x 60 mL	2 x 20 mL	160	
	23513	1 x 60 mL	1 x 20 mL	80	
PROTEIN (URINE+CSF)					
<i>Pyrogallol Red</i>					
Manual	11501	4 x 50 mL	-	200	Std Inc
	11559	2 x 250 mL	-	500	Std Inc
A15/A25 Automated System	12501	5 x 50 mL	-	250	Std Inc
BA Automated System	21512	8 x 20 mL	-	160	
	23512	4 x 20 mL	-	80	
RHEUMATOID FACTORS (RF)					
<i>Immunoturbidimetric latex assay</i>					
Manual	31922	1 x 40 mL	1 x 10 mL	50	
	31030	2 x 160 mL	2 x 40 mL	400	
A15/A25 Automated System	13922	1 x 40 mL	1 x 10 mL	50	
BA Automated System	22922	4 x 60 mL	4 x 15 mL	300	
	23922	1 x 60 mL	1 x 15 mL	75	
TOTAL BILE ACIDS					
<i>Enzymatic Cycling</i>					
Manual / A15/A25 Automated System	12551	1 x 18 mL	1 x 6 mL	24	Std Inc
BA Automated System	23551	1 x 60 mL	1 x 20 mL	80	
TOTAL IRON BINDING CAPACITY (TIBC)					
<i>Fe³⁺ / Magnesium Hydroxide Carbonate</i>					
Manual	11554	50 mL	3 g	50	
TRANSFERRIN					
<i>Immunoturbidimetric assay</i>					
Manual	31091	1 x 20 mL	-	20	
	31092	1 x 50 mL	-	50	
	31093	1 x 250 mL	-	250	
A15/A25 Automated System	13091	1 x 50 mL	-	50	
BA Automated System	22105	2 x 60 mL	2 x 15 mL	150	
	23105	1 x 60 mL	1 x 15 mL	75	
TRIGLYCERIDES					
<i>Glycerol Phosphate Oxidase/Peroxidase</i>					
Manual	11828	1 x 50 mL	-	50	Std Inc
	11528	4 x 50 mL	-	200	Std Inc
	11529	2 x 250 mL	-	500	Std Inc
A15/A25 Automated System	12528	10 x 50 mL	-	500	
BA Automated System	21528	10 x 60 mL	-	600	
	23528	4 x 60 mL	-	240	

Test	Code	RA	RB	mL/Kit	Comment
UNSATURATED IRON BINDING CAPACITY (UIBC)					
<i>Ferrozine</i>					
Manual	11835	1 x 40 mL	1 x 10 mL	50	
A15/A25 Automated System	12835	1 x 40 mL	1 x 10 mL	50	
BA Automated System	21835	2 x 60 mL	2 x 15 mL	150	
UREA/ BUN-COLOR					
<i>Urease/ Salicylate</i>					
Manual	11536	2 x 50 mL	2 x 50 mL	200	Std Inc, RA
	11537	1 x 250 mL	1 x 250 mL	500	Std Inc, RA
UREA/ BUN-UV					
<i>Urease/ Glutamate Dehydrogenase</i>					
Manual	11516	4 x 40 mL	4 x 10 mL	200	Std Inc, WR
	11517	2 x 200 mL	2 x 50 mL	500	Std Inc, WR
	11541	1 x 800 mL	1 x 200 mL	1000	Std Inc, WR
A15/A25 Automated System	12516	5 x 40 mL	5 x 10 mL	250	WR
BA Automated System	21516	8 x 60 mL	8 x 15 mL	600	
	23516	4 x 60 mL	4 x 15 mL	300	
URIC ACID					
<i>Uricase/Peroxidase</i>					
Manual	11821	1 x 50 mL	-	50	Std Inc
	11521	1 x 200 mL	-	200	Std Inc
	11522	1 x 500 mL	-	500	Std Inc
	11540	1 x 1000 mL	-	1000	Std Inc
A15/A25 Automated System	12521	10 x 50 mL	-	500	
BA Automated System	21521	10 x 60 mL	-	600	
	23521	4 x 60 mL	-	240	
ZINC					
<i>Bromo-PAPS</i>					
Manual / A15/A25 Automated System	12526	2 x 20 mL	1 x 10 mL	50	Std Inc, RA

